

# Avalo® Series Medical Cart Line

## Anesthesia Cart

The Avalo Series Standard Anesthesia Cart is designed with a synergy of function and technology to meet the unique needs of your anesthesia department. The Anesthesia cart is available in a selection of height options, drawer depths and storage configurations to streamline workflow.

Flexible divider systems, removable drawer trays, and an array of accessory options promote orderly storage and organization. A choice of (3) three innovative locking systems ensure that anesthesia medications and supplies are securely stored.



### Anesthesia Cart Features

- Precision-welded steel frame with durable, high-impact panels
- Lightweight design and profile simplify movement and steering
- A choice of innovative locking systems ensure medications and supplies are secure
- Seamless drawer trays are easy-to-clean and feature premium, self-closing drawer slides
- Wide selection of accessories and organization systems
- Expansive work area with optional slide-out work surface
- Delivered fully assembled, no tools required
- 10-year limited warranty

### Accessories & Options



Bulk Storage Bins



Sharps Container




I.V. Pole



Divided Drawer Tray



Tilt Bin Organizer



Slide-Out Shelf

Contact Capsa to inquire about additional Accessories & Options

### Cart Specifications

	Surface Height	Depth	Width
<b>Standard</b>	43.5" / 111 cm	24" / 61 cm	31" / 79 cm
<b>Intermediate</b>	39.5" / 100 cm	24" / 61 cm	31" / 79 cm
<b>Compact</b>	36" / 91 cm	24" / 61 cm	31" / 79 cm

### Drawer Specifications

	Height	Depth	Width
<b>Three Inch</b>	3" / 7.6 cm	17.5" / 45 cm	21.5" / 55 cm
<b>Six Inch</b>	6.25" / 16 cm	17.5" / 45 cm	21.5" / 55 cm
<b>Ten Inch</b>	10" / 25 cm	17.5" / 45 cm	21.5" / 55 cm

### Model Options



Standard



Intermediate



Compact

Specifications, availability & components are subject to change without notice.